



YamayBio

Ready-to-use BCA Protein Assay Kit

QUICK START GUIDE

Research use only.
Not for diagnostic procedures.

Contents and storage

Ready-to-use BCA Protein Assay Kit

BH5484 500 Assays

Storage: Store at 4°C /-20°C for 12 months. Product shipped in ice packs.

Component	Cat. No.	Volume	Storage
BCA Reagent A	BH5401	100 mL	4°C
BCA Reagent B	BH5402	3 mL	
BSA standards ① 0 µg/mL	BH5411	1 mL	-20°C
BSA standards ② 125 µg/mL	BH5413	1 mL	
BSA standards ③ 250 µg/mL	BH5416	1 mL	
BSA standards ④ 500 µg/mL	BH5419	1 mL	
BSA standards ⑤ 750 µg/mL	BH5420	1 mL	
BSA standards ⑥ 1000 µg/mL	BH5421	1 mL	
BSA standards ⑦ 1500 µg/mL	BH5422	1 mL	
BSA standards ⑧ 2000 µg/mL	BH5423	1 mL	

Introduction

The Ready-to-use BCA Protein Assay Kit is based on the BCA (bicinchoninic acid) method, enabling rapid, stable, and sensitive determination of protein concentration. In an alkaline environment, peptide bonds in proteins reduce Cu^{2+} to Cu^+ , forming a purple-colored complex with BCA reagents. This complex absorbs light strongly at 562 nm, and the absorbance is directly proportional to the protein concentration.

The kit is compatible with a wide range (**20 to 2000 µg/mL**) of samples and is tolerant to common interfering substances such as metal ions, reducing agents, chelating agents, and detergents. It also includes pre-diluted BSA standards, eliminating the need for labor-intensive dilution steps and improving accuracy and convenience.

Quick Protocol

1. Mix Reagents A and B 50:1 to create the Working Reagent.
2. **20 µL** sample or standards + **200 µL** of Working Reagent.
3. Incubate for 30 minutes at 37°C.
4. Read at 562nm.

Prepare BCA working reagent (WR)

1. Use the following formula to calculate the total volume of WR required for the assay: Total volume of WR = (# standards + # samples) × (# replicates) × (volume of WR per sample)
Example: 8 standards, 3 samples and 3 replicates : (8 standards + 3 unknowns) × (3 replicates) × (200µL) = 6.6 mL WR
2. According to the calculated amount of the required working reagent, mix the reagent A and reagent B at 50:1 volume ratio.

Note:

1. Due to the accidental error in sample addition, it is recommended to prepare extra 1-2 samples of BCA working solution to ensure a sufficient volume of WR.
2. The newly prepared BCA working solution can be stored stably for 24 hours at room temperature under sealed conditions.

Detail Procedures for BCA Protein Assay (Microplate)

1. Add **20 µL** of ready-to-use BSA standards (① to ⑧) to a 96-well microplate.
Note: Ensure BSA standards are fully dissolved and thoroughly mixed before use.
2. Add **20 µL** of each sample (in replicate) into designated wells. (Working range: 20–2000 µg/mL)
Note: It is strongly recommended to dilute samples with 1× PBS or 0.9% saline to ensure they fall within the working range.
3. Add 200 µL of WR to each well. Mix thoroughly using a plate shaker for 30 seconds. Cover the plate and incubate at 37 °C for 30 minutes.
4. Allow the plate to equilibrate to room temperature. Measure absorbance at 562 nm using a microplate reader.
5. Subtract the average 562 nm absorbance of the blank standard replicates from the readings of all other standards and unknown samples.
6. Prepare a standard curve by plotting the average blank-corrected absorbance at 562 nm for each BSA standard against its concentration (µg/mL). Use this curve to determine the protein concentration of unknown samples.

Note:

1. If using curve-fitting algorithms associated with a microplate reader, a four-parameter (quadratic) or best-fit curve provides more accurate results than a purely linear fit.
2. Data processing requires the removal of obviously incorrect values.

Compatible substance concentrations

Salts/Buffers	Compatible Concentration	Detergent	Compatible Concentration
Ammonium Sulfate	disturb	Brij35	1%
Sodium Chloride	1M	CHAPS	1%
Urea	3M	Guanidine Hydrochloride	4M
Acetate	0.2M	NP-40	1%
Glycine	1M	Siminous Glucose	1%
HEPES	0.1M	SDS	1%
MES	50mM	Triton X-100	1%
MOPS	50mM		
Sodium Citrate	<1mM	Saccharides	Compatible Concentration
PIPES	50mM	Glucose	10mM
Sodium Phosphate	0.1M	Saccharose	1M
Sodium Acetate	0.2M, pH 5.5	Reducing Agent	Compatible Concentration
TES	50mM	beta-mercaptoethanol	50µM
Tris	0.1M	DTT	1mM
Polar compound	Compatible Concentration	Other	Compatible Concentration
DMSO	5%	HCl/NaOH	0.1M
Glycerinum	10%	Lipid	disturb
Chelating Agent	Compatible Concentration		
EDTA	10mM		