

# UltraPure Mitochondria and Mitochondrial Protein Extraction Kit

## Key Features



### High Purity

The protocol uses a two-step density gradient centrifugation method, ensuring a higher purity of the isolated mitochondria.



### High Activity

The isolated mitochondria retain their intact structure and exhibit high biological activity.



### Wide Application

The mitochondria extracted using this product can be used for:

- Mitochondrial function or enzyme activity studies
- Mitochondrial DNA (mtDNA) extraction
- Mitochondrial protein analysis, such as SDS-PAGE, Western Blot (WB), Immunoprecipitation (IP), and proteomics research

Cat No.

BF7478

Size

30 Assays

## High Purity

100 mg of mouse liver tissue was processed according to the protocol to obtain the following fractions: liver homogenate, pellet (containing nuclei, unlysed cells, etc.), supernatant (cytoplasmic components), crude mitochondria, and purified mitochondria. Each fraction was diluted 1:100 in PBS and analyzed by flow cytometry.

**As shown in Fig 1, the purified mitochondrial fraction displayed a distinct, uniformly distributed cluster, indicating a high degree of mitochondrial purity.**




Fig 1. Flow cytometric analysis of different mouse liver fractions

Mitochondria extracted from 100mg mouse liver tissue using BF7478 were stained with Janus Green solution.

**Microscopic examination showed: Blue-green stained mitochondria appearing as rod-shaped or dumbbell-shaped structures, with frequent aggregation due to high mitochondrial density in liver tissue.**




Fig 2. Janus Green Staining of Mouse Liver Mitochondria (200×)

## High Activity

The isolated mitochondria were treated with MitoTracker Red, a membrane potential-dependent mitochondrial dye, and observed under a fluorescence microscope.

**As shown in Figure 3, the stained mitochondrial suspension emitted red fluorescence, indicating that the mitochondrial membrane potential was well preserved and the extracted mitochondria remained highly active.**




Fig 3. MitoTracker Staining of Mouse Liver Mitochondria




Fig 4. Western Blot Analysis of Mitochondrial Proteins

Mt: mitochondrial proteins extracted from mouse liver and HEK293 cells using BF7478.  
Total: total proteins extracted from liver tissue and HEK293 cells using RIPA lysis buffer.  
10 µg of sample was loaded per well after quantification by BCA assay.

## Wide Application

The extracted mitochondria can be used for mitochondrial function studies, mitochondrial protein analysis, and mtDNA extraction.

**As shown in Figure 4, COX IV and Hsp60 bands were much stronger in the mitochondrial fraction than in the total protein fraction, showing effective enrichment by BF7478. No Lamin A/C or Syntaxin 16 bands were detected, indicating minimal contamination from nuclei or Golgi.**