

# CERTIFICATE OF ANALYSIS

**Product Name:** Multicolor Low Range Protein Ladder 1.2-40 kDa  
**Cat.No:** CY5461-250  
**Lot.No:** A9160  
**Production Date:** 09/02/2025  
**Expiry Date:** 09/02/2027  
**Storage Temperature:** Stored at 4°C for 2 months, at -20°C for 24 months

| Parameter         | Specification   | Result   |
|-------------------|---|----------|
| Visual appearance | Blue solution   | Conforms |
| Band pattern      | 5µl of CY5461 Protein ladder was determined by performing SDS-PAGE assay. 7 clear protein bands were observed with normal band patterns and distribution. | Conforms |

Quality authorized by: *Esther Zhang*

**Issue date:** 09/05/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)