

Protein-Free Blocking Buffer(1×)

Key Features



Rapid and efficient

Complete blocking in 5 -10 minutes



Protein Free

No cross-reactivity with antibodies or target proteins, higher signal-noise ratio compared to non-fat milk and BSA. Suitable for the detection for phosphoproteins



Cat. No.

Size

BM3152 -500

500mL

Applicable to phosphorylated proteins



Higher signal-noise ratio compared to non-fat milk and BSA



Target: HSP60(Cat.No P900009), dilution at 1:5000

MW: 63kD

Lysates: Hela whole cell lysates

Secondary antibody: Goat anti-Rabbit IgG(H+L),HRP conjugated at 1:5000 dilution