

Document: RF1551-A
Release date: Sep 18, 2024

Applies to part numbers:
RF1551,RF1553,RF155

Section 1 Product and Company Identification

Product code: RF1551 Product name: AIO-Stacking Buffer A

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors, or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

The product is a mix containing no substances which, at their given concentration, are considered to be hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. Get medical attention if symptoms occur

INHALATION EXPOSURE Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. if necessary, call a physician.

DERMAL EXPOSURE Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.

EYE EXPOSURE Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.

MOST IMPORTANT SYMPTOMS AND EFFECTS,BOTH ACUTE AND DELAYED

Prolonged contact may cause redness and irritation.

INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

Treat symptomatically.

Section 5 Fire Fighting Measures

FLASH POINT N/A

AUTOIGNITION TEMP N/A

FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable	Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable	No information available
FIREFIGHTING	Not known

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES	Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.
ENVIRONMENTAL PRECAUTIONS	Avoid discharge into drains and waterways whenever possible.
METHODS AND MATERIAL FOR CONTAINMENT	Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING	Avoid inhalation, contact with eyes and skin. Always wear recommended Personal Protective Equipment.
SAFE STORAGE	Keep container tightly sealed in cool, well-ventilated area. Store in accordance with local regulations.

Section 8 Exposure Controls / Personal Protection

PERSONAL PROTECTIVE EQUIPMENT	
Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.
GENERAL HYGIENE MEASURES	Wash thoroughly after handling.
ENVIRONMENTAL EXPOSURE CONTROLS	Avoid discharge into drains and waterways whenever possible.

Section 9 Physical and Chemical Properties

Appearance Physical State	Liquid
Solvent Content	N/A

Section 10 Stability and Reactivity

STABILITY	Stable under recommended storage conditions.
HAZARDOUS DECOMPOSITION PRODUCT	
Hazardous Decomposition Products	No data available.
HAZARDOUS POLYMERIZATION	
Hazardous Polymerization	No data available.

Section 11 Toxicological Information

ROUTE OF EXPOSURE	
Skin Contact	May cause skin irritation.

Skin Absorption	May be harmful if absorbed through the skin.
Eye Contact	May cause eye irritation.
Inhalation	Material may be irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.
Ingestion	May be harmful if swallowed.

SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 Ecological Information

Toxicity	The environmental impact of this product has not been fully investigated.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	No further relevant information available.
Results of PBT and vPvB assessment	No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in according to approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG

Not regulated

ENVIRONMENTAL HAZARDS

Not Applicable

SPECIAL PRECAUTIONS FOR USER

Not Applicable

TRANSPORT IN BULK ACCORDING TO ANNEX II OF MARPOL 73/78 AND THE IBC CODE

Not Applicable

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313	This product is not regulated by SARA.
Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product does not contain any Proposition 65 chemicals.
WHMIS Hazard Class	Non-controlled

Section 16 Other Information

Date of last revision:	Sep 18,2024
------------------------	-------------

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet

Document: RF1552-A
Release date: Sep 18, 2024

Applies to part numbers:
RF1552, RF1554, RF1556

Section 1 Product and Company Identification

Product code: RF1552 Product name: AIO-Stacking Buffer B

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

GHS CLASSIFICATION

Signal Word

DANGER

Hazard pictograms



Hazard Statements

H301 - Toxic if swallowed
H411 - Toxic to aquatic life with long lasting effects
H312 + H332 - Harmful in contact with skin or if inhaled
H315 - Causes skin irritation
H319 - Causes serious eye irritation
H317 - May cause an allergic skin reaction
H372 - Causes damage to organs through prolonged or repeated exposure
H350 - May cause cancer
H340 - May cause genetic defects
H361 - Suspected of damaging fertility or the unborn child

Precautionary Statements

Prevention

P264 - Wash hands thoroughly after handling
P280 - Wear protective gloves/protective clothing/eye protection/face protection
P261 - Avoid breathing dust/fume/gas/mist/vapors/spray
P201 - Obtain special instructions before use
P260 - Do not breathe dust/fume/gas/mist/vapors/spray

Response

P308 + P313 - IF exposed or concerned: Get medical advice/attention
P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Disposal

P502 - Refer to manufacturer/supplier for information on recovery/recycling

HMIS RATING

HEALTH	2* Chronic Hazard
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

Ingredient	CAS-No.	EINECS-No.	Weight %
Acrylamide	79-06-1	201-173-7	9-30
Bis-Acrylamide	110-26-9	203-750-9	0.3-1

The product contains no other substances which, at their given concentration, are considered hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice. If swallowed, wash out mouth with water provided person is conscious. Do NOT induce vomiting. Call a physician.
INHALATION EXPOSURE	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. If necessary, call a physician.
DERMAL EXPOSURE	In case of contact, rinse with plenty of water. Remove contaminated clothing and shoes.
EYE EXPOSURE	In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.
MOST IMPORTANT SYMPTOMS AND EFFECTS, BOTH ACUTE AND DELAYED	

H301 - Toxic if swallowed H312 + H332 - Harmful in contact with skin or if inhaled H315 - Causes skin irritation H319 - Causes serious eye irritation H317 - May cause an allergic skin reaction H372 - Causes damage to organs through prolonged or repeated exposure H350 - May cause cancer H340 - May cause genetic defects H361 - Suspected of damaging fertility or the unborn child

Section 5 Fire Fighting Measures

FLASH POINT	N/A
AUTOIGNITION TEMP	N/A
FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable	Water spray, carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable	No information available
FIREFIGHTING	Not known

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES	Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.
ENVIRONMENTAL PRECAUTIONS	Avoid discharge into drains and waterways whenever possible.
METHODS AND MATERIAL FOR CONTAINMENT	Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING	Avoid inhalation, contact with eyes and skin. Always wear recommended Personal Protective Equipment.
---------------	--

SAFE STORAGE

Keep container tightly sealed in cool, well-ventilated area. Store in accordance with local regulations.

Section 8 Exposure Controls / Personal Protection**PERSONAL PROTECTIVE EQUIPMENT**

Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.

GENERAL HYGIENE MEASURES

Wash thoroughly after handling.

ENVIRONMENTAL EXPOSURE CONTROLS

Avoid discharge into drains and waterways whenever possible.

Section 9 Physical and Chemical Properties

Appearance Physical State	Liquid
Solvent Content	N/A

Section 10 Stability and Reactivity

STABILITY Stable under recommended storage conditions.

HAZARDOUS DECOMPOSITION PRODUCT

Hazardous Decomposition Products No data available.

HAZARDOUS POLYMERIZATION

Hazardous Polymerization No data available.

Section 11 Toxicological Information**INFORMATION ON TOXICOLOGICAL EFFECTS**

Ingredient	Oral LD50	Dermal LD50	Inhalation LC50
Acrylamide	= 177 mg/kg	= 1148 mg/kg	No data available
Bis-Acrylamide	= 390 mg/kg (Rat)	No data available	No data available

PRINCIPAL ROUTES OF EXPOSURE

Skin corrosion/irritation	Skin irritation Serious eye damage/eye irritation
Respiratory or skin sensitization	May cause sensitization by skin contact
Specific target organ toxicity (STOT)–repeated exposure	Target organ(s): Blood
Carcinogenicity	Contains a known or suspected carcinogen
Germ cell mutagenicity	Product is or contains a chemical which may cause mutations in the germ cells of humans.
Reproductive toxicity	May cause adverse reproductive effects - such as birth defect, miscarriages, or infertility

Section 12 Ecological Information**TOXICITY**

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Acrylamide	No data available	Daphnia magna EC50=98 mg/L (48 h)	No data available	No data available	logPow-0.9
Bis-Acrylamide	No data available	No data available	No data available	No data available	logPow-0.08

Persistence and degradability No further relevant information available.

Bioaccumulative potential No further relevant information available.

General notes No further relevant information available.

Results of PBT and vPvB assessment No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG

Not regulated

ENVIRONMENTAL HAZARDS

Not Applicable

SPECIAL PRECAUTIONS FOR USER

Not Applicable

TRANSPORT IN BULK ACCORDING TO ANNEX II OF MARPOL 73/78 AND THE IBC CODE

Not Applicable

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313

This product contains the following toxic chemical(s) subject to the notification requirements of section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986. This law requires certain manufacturers to report on annual emissions of specified chemicals and chemical categories. Please note that if you repackage, or otherwise redistribute, this product to industrial customers, a notice similar to this one should be sent to those customers.

Ingredient	CAS #	Weight %	SARA 313 - Threshold Values%
Acrylamide	7758-99-8	9-30	0.1

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product contains the following HAPs

Ingredient	CAS #	Weight %	HAPS data
Acrylamide	7758-99-8	9-30	Present

US STATE REGULATIONS

California Proposition 65

This product contains the following Proposition 65 chemicals:

Ingredient	CAS #	Weight %	Category
Acrylamide	7758-99-8	9-30	Carcinogen Developmental Male Reproductive

WHMIS Hazard Class

D1B - Toxic materials

Section 16 Other Information

Date of last revision: Sep 18,2024

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet

Document: RF1553-A
Release date: Sep 18, 2024

Applies to part numbers:
RF1551,RF1553,RF1555

Section 1 Product and Company Identification

Product code: RF1553 Product name: AIO-Resolving Buffer A (10%)

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors, or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

The product is a mix containing no substances which, at their given concentration, are considered to be hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. Get medical attention if symptoms occur

INHALATION EXPOSURE Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. if necessary, call a physician.

DERMAL EXPOSURE Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.

EYE EXPOSURE Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.

MOST IMPORTANT SYMPTOMS AND EFFECTS,BOTH ACUTE AND DELAYED

Prolonged contact may cause redness and irritation.

INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

Treat symptomatically.

Section 5 Fire Fighting Measures

FLASH POINT N/A

AUTOIGNITION TEMP N/A

FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable	Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable	No information available
FIREFIGHTING	Not known

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES	Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.
ENVIRONMENTAL PRECAUTIONS	Avoid discharge into drains and waterways whenever possible.
METHODS AND MATERIAL FOR CONTAINMENT	Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING	Avoid inhalation, contact with eyes and skin. Always wear recommended Personal Protective Equipment.
SAFE STORAGE	Keep container tightly sealed in cool, well-ventilated area. K. Store in accordance with local regulations.

Section 8 Exposure Controls / Personal Protection

PERSONAL PROTECTIVE EQUIPMENT	
Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.
GENERAL HYGIENE MEASURES	Wash thoroughly after handling.
ENVIRONMENTAL EXPOSURE CONTROLS	Avoid discharge into drains and waterways whenever possible.

Section 9 Physical and Chemical Properties

Appearance Physical State	Liquid
Solvent Content	N/A

Section 10 Stability and Reactivity

STABILITY	Stable under recommended storage conditions.
HAZARDOUS DECOMPOSITION PRODUCT	
Hazardous Decomposition Products	No data available.
HAZARDOUS POLYMERIZATION	
Hazardous Polymerization	No data available.

Section 11 Toxicological Information

ROUTE OF EXPOSURE	
Skin Contact	May cause skin irritation.

Skin Absorption	May be harmful if absorbed through the skin.
Eye Contact	May cause eye irritation.
Inhalation	Material may be irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.
Ingestion	May be harmful if swallowed.

SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 Ecological Information

Toxicity	The environmental impact of this product has not been fully investigated.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	No further relevant information available.
Results of PBT and vPvB assessment	No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in according to approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG

Not regulated

ENVIRONMENTAL HAZARDS

Not Applicable

SPECIAL PRECAUTIONS FOR USER

Not Applicable

TRANSPORT IN BULK ACCORDING TO ANNEX II OF MARPOL 73/78 AND THE IBC CODE

Not Applicable

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313	This product is not regulated by SARA.
Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product does not contain any Proposition 65 chemicals.
WHMIS Hazard Class	Non-controlled

Section 16 Other Information

Date of last revision:	Sep 18,2024
------------------------	-------------

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet

Document: RF1554-A
Release date: Sep 18, 2024

Applies to part numbers:
RF1552, RF1554, RF1556

Section 1 Product and Company Identification

Product code: RF1554 Product name: AIO-Resolving Buffer B (10%)

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

GHS CLASSIFICATION

Signal Word

DANGER

Hazard pictograms



Hazard Statements

H301 - Toxic if swallowed
H411 - Toxic to aquatic life with long lasting effects
H312 + H332 - Harmful in contact with skin or if inhaled
H315 - Causes skin irritation
H319 - Causes serious eye irritation
H317 - May cause an allergic skin reaction
H372 - Causes damage to organs through prolonged or repeated exposure
H350 - May cause cancer
H340 - May cause genetic defects
H361 - Suspected of damaging fertility or the unborn child

Precautionary Statements

Prevention

P264 - Wash hands thoroughly after handling
P280 - Wear protective gloves/protective clothing/eye protection/face protection
P261 - Avoid breathing dust/fume/gas/mist/vapors/spray
P201 - Obtain special instructions before use
P260 - Do not breathe dust/fume/gas/mist/vapors/spray

Response

P308 + P313 - IF exposed or concerned: Get medical advice/attention
P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Disposal

P502 - Refer to manufacturer/supplier for information on recovery/recycling

HMIS RATING

HEALTH	2* Chronic Hazard
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

Ingredient	CAS-No.	EINECS-No.	Weight %
Acrylamide	79-06-1	201-173-7	9-30
Bis-Acrylamide	110-26-9	203-750-9	0.3-1

The product contains no other substances which, at their given concentration, are considered hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice. If swallowed, wash out mouth with water provided person is conscious. Do NOT induce vomiting. Call a physician.
INHALATION EXPOSURE	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. If necessary, call a physician.
DERMAL EXPOSURE	In case of contact, rinse with plenty of water. Remove contaminated clothing and shoes.
EYE EXPOSURE	In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.
MOST IMPORTANT SYMPTOMS AND EFFECTS, BOTH ACUTE AND DELAYED	

H301 - Toxic if swallowed H312 + H332 - Harmful in contact with skin or if inhaled H315 - Causes skin irritation H319 - Causes serious eye irritation H317 - May cause an allergic skin reaction H372 - Causes damage to organs through prolonged or repeated exposure H350 - May cause cancer H340 - May cause genetic defects H361 - Suspected of damaging fertility or the unborn child

Section 5 Fire Fighting Measures

FLASH POINT	N/A
AUTOIGNITION TEMP	N/A
FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable	Water spray, carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable	No information available
FIREFIGHTING	Not known

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES	Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.
ENVIRONMENTAL PRECAUTIONS	Avoid discharge into drains and waterways whenever possible.
METHODS AND MATERIAL FOR CONTAINMENT	Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING	Avoid inhalation, contact with eyes and skin. Always wear recommended Personal Protective Equipment.
---------------	--

SAFE STORAGE

Keep container tightly sealed in cool, well-ventilated area. Store in accordance with local regulations.

Section 8 Exposure Controls / Personal Protection**PERSONAL PROTECTIVE EQUIPMENT**

Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.

GENERAL HYGIENE MEASURES

Wash thoroughly after handling.

ENVIRONMENTAL EXPOSURE CONTROLS

Avoid discharge into drains and waterways whenever possible.

Section 9 Physical and Chemical Properties

Appearance Physical State	Liquid
Solvent Content	N/A

Section 10 Stability and Reactivity

STABILITY Stable under recommended storage conditions.

HAZARDOUS DECOMPOSITION PRODUCT

Hazardous Decomposition Products No data available.

HAZARDOUS POLYMERIZATION

Hazardous Polymerization No data available.

Section 11 Toxicological Information**INFORMATION ON TOXICOLOGICAL EFFECTS**

Ingredient	Oral LD50	Dermal LD50	Inhalation LC50
Acrylamide	= 177 mg/kg	= 1148 mg/kg	No data available
Bis-Acrylamide	= 390 mg/kg (Rat)	No data available	No data available

PRINCIPAL ROUTES OF EXPOSURE

Skin corrosion/irritation	Skin irritation Serious eye damage/eye irritation
Respiratory or skin sensitization	May cause sensitization by skin contact
Specific target organ toxicity (STOT)–repeated exposure	Target organ(s): Blood
Carcinogenicity	Contains a known or suspected carcinogen
Germ cell mutagenicity	Product is or contains a chemical which may cause mutations in the germ cells of humans.
Reproductive toxicity	May cause adverse reproductive effects - such as birth defect, miscarriages, or infertility

Section 12 Ecological Information**TOXICITY**

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Acrylamide	No data available	Daphnia magna EC50=98 mg/L (48 h)	No data available	No data available	logPow-0.9
Bis-Acrylamide	No data available	No data available	No data available	No data available	logPow-0.08

Persistence and degradability No further relevant information available.

Bioaccumulative potential No further relevant information available.

General notes No further relevant information available.

Results of PBT and vPvB assessment No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG

Not regulated

ENVIRONMENTAL HAZARDS

Not Applicable

SPECIAL PRECAUTIONS FOR USER

Not Applicable

TRANSPORT IN BULK ACCORDING TO ANNEX II OF MARPOL 73/78 AND THE IBC CODE

Not Applicable

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313

This product contains the following toxic chemical(s) subject to the notification requirements of section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986. This law requires certain manufacturers to report on annual emissions of specified chemicals and chemical categories. Please note that if you repackage, or otherwise redistribute, this product to industrial customers, a notice similar to this one should be sent to those customers.

Ingredient	CAS #	Weight %	SARA 313 - Threshold Values%
Acrylamide	7758-99-8	9-30	0.1

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product contains the following HAPs

Ingredient	CAS #	Weight %	HAPS data
Acrylamide	7758-99-8	9-30	Present

US STATE REGULATIONS

California Proposition 65

This product contains the following Proposition 65 chemicals:

Ingredient	CAS #	Weight %	Category
------------	-------	----------	----------

Acrylamide	7758-99-8	9-30	Carcinogen Developmental Male Reproductive
------------	-----------	------	--

WHMIS Hazard Class

D1B - Toxic materials

Section 16 Other Information

Date of last revision: Sep 18,2024

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet

Document: RF1557-A
Release date: Aug 18, 2024

Applies to part numbers:
RF1557,GB4353-5,GB4353-15

Section 1 Product and Company Identification

Product code: RF1557 Product name: Tricine Protein Sample Loading Buffer (Denaturing, Reducing, 2×)

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors, or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

Ingredient	CAS-No.	EINECS-No.	Weight %
Glycerol	56-81-5	200-289-5	40-70
Lauryl alcohol sulfate	2044-56-6	218-055-6	2-5

Section 4 First Aid Measures

ORAL EXPOSURE If swallowed, wash out mouth with water provided person is conscious. Call a physician.
INHALATION EXPOSURE If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
DERMAL EXPOSURE In case of contact, immediately wash skin with soap and copious amounts of water.
EYE EXPOSURE In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

Section 5 Fire Fighting Measures

FLASH POINT N/A
AUTOIGNITION TEMP N/A
FLAMMABILITY N/A
EXTINGUISHING MEDIA
Suitable Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable No information available
FIREFIGHTING Standard procedure for chemical fires

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES Ensure adequate ventilation.
Always wear recommended Personal Protective Equipment.
Use personal protection equipment .

ENVIRONMENTAL PRECAUTIONS No special environmental precautions required.

METHODS AND MATERIAL FOR CONTAINMENT Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING Use personal protective equipment as required. No special handling advices are necessary.

SAFE STORAGE Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers.

Section 8 Exposure Controls / Personal Protection

PERSONAL PROTECTIVE EQUIPMENT

Respiratory In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

Eye Tight sealing safety goggles.

Skin and Body Wear suitable protective clothing.

GENERAL HYGIENE MEASURES

Wash thoroughly after handling.

ENVIRONMENTAL EXPOSURE CONTROLS

No special environmental precautions required.

Section 9 Physical and Chemical Properties

Appearance Physical State Liquid

Solvent Content N/A

Section 10 Stability and Reactivity

STABILITY Stable under normal conditions.

HAZARDOUS DECOMPOSITION PRODUCT

Hazardous Decomposition Products No dangerous decomposition products known

HAZARDOUS POLYMERIZATION

Hazardous Polymerization No data available.

Section 11 Toxicological Information

Information on toxicological effects

Ingredient	Oral LD50	Dermal LD50	Inhalation LC50
Glycerol	= 12600 mg/kg Oral	> 10 g/kg(Rabbit)	>2.75mg/L(Rat)
Lauryl alcohol sulfate	No data available.	No data available	No data available

ROUTE OF EXPOSURE

Skin Contact May cause skin irritation.

Skin Absorption May be harmful if absorbed through the skin.

Eye Contact May cause eye irritation.

Inhalation	Material may be irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.
Ingestion	May be harmful if swallowed.

SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 Ecological Information

Toxicity	The environmental impact of this product has not been fully investigated.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water. Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment	Not Applicable.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

Observe all federal, state, and local environmental regulations.

Within the present knowledge of the supplier, this product is not regarded as hazardous waste.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG	Not regulated in the meaning of transport regulations.
Environmental hazards	Not Applicable
Special precautions for user	Not Applicable
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	Not Applicable

Section 15 Regulatory Information

Ingredient	US TSCA
Glycerol 56-81-5 (40-70)	Listed
Lauryl alcohol sulfate, lithium salt 2044-56-6 (<10)	Listed

US FEDERAL REGULATIONS

SARA 313	This product is not regulated by SARA.
Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product does not contain any Proposition 65 chemicals.
WHMIS Hazard Class	Non-controlled

Section 16 Other Information

Date of last revision:	Aug 18,2024
------------------------	-------------

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet



Document: MS8136-A
Release date: Aug 18, 2024

Applies to part numbers:
MS8136

Section 1 Product and Company Identification

Product code: MS8136 Product name: Tris/Tricine/SDS Buffer Instant Granules
Supplier: YAMAYBIO INC
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA
Recommended Use: For research use only.
Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors, or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

Chemical characterization: Mixtures

Ingredient	CAS No.	EINECS-No	Weight %
sodium dodecyl sulfate	151-21-3	205-788-1	3-6

There are no other ingredients present which within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous according to the OSHA Hazard Communications Standard (29 CFR 1910.1200) and require reporting in this section.

Section 4 First Aid Measures

ORAL EXPOSURE If swallowed, wash out mouth with water provided person is conscious. Call a physician.
INHALATION EXPOSURE If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
DERMAL EXPOSURE In case of contact, wash skin with water.
EYE EXPOSURE In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers.

Section 5 Fire Fighting Measures

FLASH POINT N/A
AUTOIGNITION TEMP N/A
FLAMMABILITY N/A
EXTINGUISHING MEDIA
Suitable Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable No information available

FIREFIGHTING

Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

Section 6 Accidental Release Measures**PERSONAL PRECAUTIONS,
PROTECTIVE EQUIPMENT AND
EMERGENCY PROCEDURES**

Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment.

ENVIRONMENTAL PRECAUTIONS

No special environmental precautions required.

**METHODS AND MATERIAL FOR
CONTAINMENT**

Soak up with inert absorbent material.

Section 7 Handling and Storage**SAFE HANDLING**

Use personal protective equipment as required. No special handling advices are necessary.

SAFE STORAGE

Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers. Store in accordance with local regulations.

Section 8 Exposure Controls / Personal Protection**PERSONAL PROTECTIVE EQUIPMENT****Respiratory**

In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand

Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

Eye

Tight sealing safety goggles.

Skin and Body

Wear suitable protective clothing.

GENERAL HYGIENE MEASURES

Wash thoroughly after handling.

**ENVIRONMENTAL EXPOSURE
CONTROLS**

No special environmental precautions required.

Section 9 Physical and Chemical Properties**Appearance Physical State**

White or off-white powder

Solvent Content

N/A

Section 10 Stability and Reactivity**STABILITY**

Stable under normal conditions.

Materials to Avoid

Strong oxidizing agents, metals, acids.

HAZARDOUS DECOMPOSITION PRODUCT**Hazardous Decomposition
Products**

No dangerous decomposition products known

HAZARDOUS POLYMERIZATION**Hazardous Polymerization**

No data available.

Section 11 Toxicological Information**INFORMATION ON TOXICOLOGICAL EFFECTS**

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
sodium dodecyl sulfate	= 1288 mg/kg (Rat)	= 200 mg/kg(Rabbit)	>3900mg/m3(Rat)

ROUTE OF EXPOSURE**Skin Contact**

May cause skin irritation.

Skin Absorption	May be harmful if absorbed through the skin.
Eye Contact	May cause eye irritation.
Inhalation	Material may be irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.
Ingestion	May be harmful if swallowed.

SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 Ecological Information

Toxicity	The environmental impact of this product has not been fully investigated.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water. Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment	No information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in according to approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG	Not regulated in the meaning of transport regulations.
Environmental hazards	N/A
Special precautions for user	N/A
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	N/A

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313	This product is not regulated by SARA.
Clean Air Act, Section 112	
Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product does not contain any Proposition 65 chemicals.
WHMIS Hazard Class	Non-controlled

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

Section 16 Other Information

Date of last revision: Aug 18,2024

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet

Document: BM31531S-A
Release date: Sep 18, 2024

Applies to part numbers:
BM3151, BM3151-S

Section 1 Product and Company Identification

Product code: BM3151-S Product name: Rapid Coomassie Blue Staining (Destaining-free)

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/ Information on Ingredients

Ingredient	CAS No.	EINECS-No	Weight %
Ethyl alcohol	64-17-5	200-578-6	3-5

The product contains no other substances which, at their given concentration, are considered hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice. If swallowed, wash out mouth with water provided person is conscious. Call a physician.

INHALATION EXPOSURE Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. If necessary, call a physician.

DERMAL EXPOSURE In case of contact, rinse with plenty of water.

EYE EXPOSURE In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

Section 5 Fire Fighting Measures

FLASH POINT N/A

AUTOIGNITION TEMP N/A

FLAMMABILITY N/A

EXTINGUISHING MEDIA

Suitable Water spray, carbon dioxide, dry chemical powder, or appropriate foam.

Not suitable No information available.

FIREFIGHTING

Protective Equipment	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.
Specific Hazard(s)	Emits toxic fumes under fire conditions.

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.

ENVIRONMENTAL PRECAUTIONS Try to prevent further leakage or spillage. Keep the product away from drains or water courses.

METHODS AND MATERIAL FOR CONTAINMENT Absorb solutions with finely-powdered liquid-binding material (diatomite, universal binders); Decontaminate surfaces and equipment by scrubbing with alcohol; Dispose of contaminated material according to Section 13.

Section 7 Handling and Storage

SAFE HANDLING Avoid inhalation, contact with eyes and skin. Avoid dust and aerosol formation. Use only in areas with appropriate exhaust ventilation.

SAFE STORAGE Keep container tightly sealed in cool, well-ventilated area. Keep away from direct sunlight and sources of ignition.

Section 8 Exposure Controls / Personal Protection

PERSONAL PROTECTIVE EQUIPMENT

Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.

GENERAL HYGIENE MEASURES Wash thoroughly after handling.

ENVIRONMENTAL EXPOSURE CONTROLS Keep the product away from drains, water courses or the soil. Clean spillages in a safe way as soon as possible.

Section 9 Physical and Chemical Properties

APPEARANCE PHYSICAL STATE	blue liquid
SOLVENT CONTENT	N/A

Section 10 Stability and Reactivity

STABILITY Stable under recommended storage conditions.

Materials to Avoid Strong acids/alkalis, strong oxidising/reducing agents.

HAZARDOUS DECOMPOSITION PRODUCT

Hazardous Decomposition Products	Under fire conditions, may decompose and emit toxic fumes
----------------------------------	---

HAZARDOUS POLYMERIZATION

Hazardous Polymerization	No data available.
--------------------------	--------------------

Section 11 Toxicological Information

INFORMATION ON TOXICOLOGICAL EFFECTS

Ingredient	Oral LD50	Dermal LD50	Inhalation LC50
-------------------	------------------	--------------------	------------------------

Ethyl alcohol	= 7060 mg/kg (Rat)	No data available	=124.7mg/L(Rat)
---------------	--------------------	-------------------	-----------------

ROUTE OF EXPOSURE

Skin Contact	May cause skin irritation.
Eye Contact	May cause eye irritation.
Inhalation	May be harmful if inhaled. Causes respiratory tract irritation.
Ingestion	Harmful if swallowed.

SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Section 12 Ecological Information

Toxicity	No further relevant information available.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	No further relevant information available.
Results of PBT and vPvB assessment	No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. Observe all federal, state, and local environmental regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG	Not regulated in the meaning of transport regulations.
Environmental hazards	N/A
Special precautions for user	N/A
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	N/A

Section 15 Regulatory Information

Ingredient	US TSCA
Ethyl alcohol 64-17-5 (3-5)	Listed

US FEDERAL REGULATIONS

SARA 313	This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.
Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product contains the following Proposition 65 chemicals.
---------------------------	---

Chemical Name	CAS-No	Weight %	Category
Ethyl alcohol	64-17-5	3-5	Carcinogen Developmental

WHMIS Hazard Class	Non-controlled
--------------------	----------------

Section 16 Other Information

Date of last revision: Sep 18,2024

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet