

Document: BF7471-A
Release date: Sep 18, 2024

Applies to part numbers:
BF7471

Section 1 Product and Company Identification

Product code: BF7471 Product name: RIPA Lysis Buffer (Strong)

Supplier: YAMAYBIO INC

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

Section 2 Hazards Identification

EMERGENCY OVERVIEW:

Harmful if swallowed. May cause irritation. Avoid breathing vapors or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING

HEALTH	0
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

Section 3 Composition/Information on Ingredients

Ingredient	CAS-No.	EINECS-No.	Weight %
Triton X-100	9002-93-1	618-344-0	0.1-1
Sodium Deoxycholate	302-95-4	206-132-7	0.1-1

The product contains no other substances which, at their given concentration, are considered hazardous to health.

Section 4 First Aid Measures

ORAL EXPOSURE	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice. If swallowed, wash out mouth with water provided person is conscious. Call a physician.
INHALATION EXPOSURE	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. If necessary, call a physician.
DERMAL EXPOSURE	In case of contact, rinse with plenty of water.
EYE EXPOSURE	In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

Section 5 Fire Fighting Measures

FLASH POINT	N/A
AUTOIGNITION TEMP	N/A
FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable	Water spray, Carbon dioxide, dry chemical powder, or appropriate foam.
Not suitable	No information available

Section 6 Accidental Release Measures

PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES	Ensure adequate ventilation. Always wear recommended Personal Protective Equipment. Use personal protection equipment. See Section 8 for more detail.
ENVIRONMENTAL PRECAUTIONS	No special environmental precautions required.
METHODS AND MATERIAL FOR CONTAINMENT	Soak up with inert absorbent material.

Section 7 Handling and Storage

SAFE HANDLING	Use personal protective equipment as required. No special handling advices are necessary.
SAFE STORAGE	Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers.

Section 8 Exposure Controls / Personal Protection

PERSONAL PROTECTIVE EQUIPMENT	
Respiratory	In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.
Hand	Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.
Eye	Tight sealing safety goggles.
Skin and Body	Wear suitable protective clothing.
GENERAL HYGIENE MEASURES	Wash thoroughly after handling.
ENVIRONMENTAL EXPOSURE CONTROLS	No special environmental precautions required.

Section 9 Physical and Chemical Properties

APPEARANCE PHYSICAL STATE	Liquid, water-based, pH 6.5 to 8
SOLVENT CONTENT	N/A

Section 10 Stability and Reactivity

STABILITY	Stable under normal conditions.
Materials to Avoid	No dangerous reaction known under conditions of normal use.
HAZARDOUS DECOMPOSITION PRODUCT	
Hazardous Decomposition Products	No dangerous decomposition products known.
HAZARDOUS POLYMERIZATION	
Hazardous Polymerization	No data available.

Section 11 Toxicological Information

ROUTE OF EXPOSURE	
Skin Contact	May cause skin irritation.
Skin Absorption	May be harmful if absorbed through the skin.
Eye Contact	May cause eye irritation.
Inhalation	Material may be irritating to mucous membranes and upper respiratory tract. May be harmful if inhaled.
Ingestion	May be harmful if swallowed.
SIGNS AND SYMPTOMS OF EXPOSURE	
To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.	

Section 12 Ecological Information

Toxicity	The environmental impact of this product has not been fully investigated.
Persistence and degradability	No further relevant information available.
Bioaccumulative potential	No further relevant information available.
General notes	No further relevant information available.
Results of PBT and vPvB assessment	No further relevant information available.

Section 13 Disposal Considerations

APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. Observe all federal, state, and local environmental regulations.

Section 14 Transport Information

IATA / ADR / DOT-US / IMDG	Not regulated in the meaning of transport regulations.
Environmental hazards	N/A
Special precautions for user	N/A
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	N/A

Section 15 Regulatory Information

US FEDERAL REGULATIONS

SARA 313	This product is not regulated by SARA.
Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)	This product does not contains HAPs.

US STATE REGULATIONS

California Proposition 65	This product does not contain any Proposition 65 chemicals.
WHMIS Hazard Class	Non-controlled

Section 16 Other Information

Date of last revision: Sep 18,2024

The information provided herein has been obtained through diligent research and analysis. The recommendations are based upon the prudent application of professional judgment. This information is not intended to be exhaustive and should be used as a general guideline. All materials and mixtures may present potential hazards, and should be handled with care. As the Company cannot control the specific methods, quantities, or conditions of use, it shall not be liable for any damages or losses resulting from the handling or contact with the product as described herein.

THE INFORMATION CONTAINED IN THIS SAFETY DATA SHEET DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

End of Safety Data Sheet