

CERTIFICATE OF ANALYSIS

Product Name: One-Step Stain-free PAGE Gel Rapid Preparation Kit (7.5%)
Cat.No: RF1522
Lot.No: A21100
Production Date: 02/01/2025
Expiry Date: 08/01/2026
Storage Temperature: Store at 4°C, the enhanced catalyst is stored at -20°C.

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
pH	8.6-9.0	Conforms
Gelation time	10-30min	Conforms
Migration	the band of 35kd of marker was seperated	Conforms
Gel visualization	5µg, 10µg, 20µg and 30µg of total proteins were loaded, then the gels were excited at 302 nm UV light, protein bands can be visualized in 5mins	Conforms

Quality authorized by:

Esther Zhang

Issue date: 02/03/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com