



YamayBio

TBS/Tween Buffer Instant Granules

INSTRUMENT MANUAL

Research use only.
Not for diagnostic procedures.

Contents and storage

TBS/Tween Buffer Instant Granules

MS8133 10 sticks (1 L×10)

Storage: Store at room temperature for 36 months. Product shipped at room temperature.

Introduction

This product is provided in convenient white or off-white instant granule sticks. Each stick contains a blend of sodium chloride, potassium chloride, Tris, and Tween-20. When dissolved, one stick yields 1 liter of 1× TBST buffer with a final composition of 137 mM sodium chloride, 2.7 mM potassium chloride, 25 mM Tris, and 0.05% Tween-20.

Quick Start Protocol

Dissolve one stick in distilled water and adjust the volume to **1 L** to prepare 1× TBS/Tween Buffer.