



YamayBio

Tris/Glycine/SDS Buffer Instant Granules

INSTRUMENT MANUAL

Research use only.
Not for diagnostic procedures.

Contents and storage

Tris/Glycine/SDS Buffer Instant Granules

MS8131

10 sticks (1 L×10)

Storage: Store at room temperature for 36 months. Product shipped at room temperature.

Introduction

This product is provided as white or pale yellow instant granules. Its main components are Tris, glycine, and SDS. Each stick yields 1 L of 1× Tris/glycine/SDS electrophoresis buffer, containing 25 mM Tris, 192 mM glycine, and 0.1% SDS.

Quick Start Protocol

Dissolve one stick in distilled water and adjust the volume to **1 L** to prepare 1× Tris/glycine/SDS electrophoresis buffer.