



YamayBio

Tris/MOPS/SDS Buffer Instant Granules

INSTRUMENT MANUAL

Research use only.
Not for diagnostic procedures.

Contents and storage

Tris/MOPS/SDS Buffer Instant Granules

MS8132-10 10 sticks (500mL×10)

MS8132-20 20 sticks (500mL×20)

Storage: Store at room temperature for 36 months. Product shipped at room temperature.

Introduction

This product is provided as white or pale yellow instant granules. The main components are Tris, MOPS, EDTA, and SDS. Each stick yields 500 mL of 1× Tris/MOPS/SDS electrophoresis buffer, containing 50 mM Tris, 50 mM MOPS, 1 mM EDTA, and 0.1% SDS.

Quick Start Protocol

Dissolve one stick in distilled water and adjust the volume to **500mL** to prepare 1× Tris/MOPS/SDS electrophoresis buffer.