

CERTIFICATE OF ANALYSIS

Product Name: Tris/MOPS/SDS Buffer Instant Granules
Cat.No: MS8132-20
Lot.No: 28C1800
Production Date: 01/22/2025
Expiry Date: 01/22/2028
Storage Temperature: Stored at room temperature

Parameter	Specification	Result
Visual appearance	white or light yellow granules	Conforms
Solubility	200s(500rpm)	Conforms
Weight	8.75-8.95g/pack	Conforms
pH	7.6-7.8	Conforms
Conductivity	3000-4250 μ S/cm	Conforms
Dnase	None	Conforms
Protease	None	Conforms
Function test	The detected bands have high resolution and clarity	Conforms

Quality authorized by:

Esther Zhang

Issue date: 01/22/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com