

## CERTIFICATE OF ANALYSIS

**Product Name:** Tris/MOPS/SDS Buffer Instant Granules  
**Cat.No:** MS8132-10  
**Lot.No:** 28B9700  
**Production Date:** 11/20/2024  
**Expiry Date:** 11/20/2027  
**Storage Temperature:** Stored at room temperature

Parameter	Specification	Result
Visual appearance	white or light yellow granules	Conforms
Solubility	200s(500rpm)	Conforms
Weight	8.75-8.95g/pack	Conforms
pH	7.6-7.8	Conforms
Conductivity	3000-4250 $\mu$ S/cm	Conforms
Dnase	None	Conforms
Protease	None	Conforms
Function test	The detected bands have high resolution and clarity	Conforms

Quality authorized by:

*Esther Zhang*

**Issue date:** 11/20/2024

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)