

CERTIFICATE OF ANALYSIS

Product Name: TBS/Tween Buffer Instant Granules
Cat.No: MS8133-10
Lot.No: 28C5850
Production Date: 12/09/2024
Expiry Date: 12/09/2027
Storage Temperature: Stored at room temperature

Parameter	Specification	Result
Visual appearance	white or off-white granules	Conforms
Solubility	100s(500rpm)	Conforms
Weight	11.89-13.15g/pack	Conforms
pH	7.25-7.55	Conforms
Conductivity	13.0-16.5 μ S/cm	Conforms
Dnase	None	Conforms
Protease	None	Conforms

Quality authorized by:

Esther Zhang

Issue date: 12/09/2024

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com