

# CERTIFICATE OF ANALYSIS

**Product Name:** Rapid Coomassie Blue Staining (Destaining-free)  
**Cat.No:** BM3151-500  
**Lot.No:** A31020  
**Production Date:** 03/05/2025  
**Expiry Date:** 03/05/2026  
**Storage Temperature:** Store at 4°C away from light

Parameter	Specification	Result
Visual appearance	Dark blue, free of particulates	Conforms
Function test	A series of bovine albumin (10ng, 20ng, 40ng, 80ng) were separated by electrophoresis on a 12.5% Tris-glycine SDS -PAGE, 1.0mm, 10wells and then stained with BM3151 Rapid Coomassie Blue Staining Solution	10ng BSA was detected

Quality authorized by:

*Esther Zhang*

**Issue date:** 03/05/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)