

CERTIFICATE OF ANALYSIS

Product Name: Tricine SDS Running Buffer Instant Granules
Cat.No: MS8136
Lot.No: 28C5876
Production Date: 12/02/2025
Expiry Date: 12/02/2028
Storage Temperature: Stored at room temperature

Parameter	Specification	Result
Visual appearance	White or off-white granules	Conforms
Solubility	100s(500rpm)	Conforms
Weight	15.50-16.05g/pack	Conforms
pH	8.0-8.3	Conforms
Conductivity	400.0-500.0 μ S/cm	Conforms
Dnase	None	Conforms
Protease	None	Conforms

Quality authorized by:

Esther Zhang

Issue date: 12/02/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com