



YamayBio

Tricine SDS Running Buffer Instant Granules

INSTRUMENT MANUAL

Research use only.
Not for diagnostic procedures.

Contents and storage

Tricine SDS Running Buffer Instant Granules

MS8136 10 sticks (500mL×10)

Storage: Store at room temperature for 36 months. Product shipped at room temperature.

Introduction

This product is provided as white or pale yellow instant granules. Each stick yields 500 mL of 1× Tricine SDS electrophoresis buffer.

When using this product for Tricine-SDS-PAGE, there is no need to prepare separate anode and cathode buffers.

Quick Start Protocol

Dissolve one stick in distilled water and adjust the volume to **500mL** to prepare 1× Tricine SDS electrophoresis buffer.