

# CERTIFICATE OF ANALYSIS

**Product Name:** ECL Western Blotting Sensitivity Substrate (Pico)  
**Cat.No:** VQ4727-100  
**Lot.No:** 55AC1020  
**Production Date:** 12/18/2025  
**Expiry Date:** 12/18/2026  
**Storage Temperature:** Store at 4°C away from light

Parameter	Specification	Result
Visual appearance	Transparent, free of particulates	Conforms
Function test	Gradient-diluted HRP was added to a 96-well black plate, followed by the addition of VQ4727 working solution, resulting in clearly developed gradient results.	Conforms

Quality authorized by:



**Issue date:** 12/18/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)