

# CERTIFICATE OF ANALYSIS

**Product Name:** Unstained Protein Ladder  
**Cat.No:** CY5462-500  
**Lot.No:** 75B1115  
**Production Date:** 01/21/2026  
**Expiry Date:** 01/21/2028  
**Storage Temperature:** Stored at 4°C for 3 months, at -20°C for 2 years

Parameter	Specification	Result
Visual appearance	Blue solution	Conforms
Band pattern	5µl of CY5462 Unstained Protein ladder was determined by performing SDS-PAGE assay. After staining with Coomassie Brilliant Blue, 10 clear protein bands were observed with normal band patterns and distribution.	Conforms

Quality authorized by:

*Esther Zhang*

**Issue date:** 01/21/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)