



YamayBio

Tris/MOPS/SDS Buffer Instant Granules

USER MANUAL

For research use only.
Not for use in diagnostic procedures.



Contents and Storage

Tris/MOPS/SDS Buffer Instant Granules

MS8132-10 10 sticks (500 mL×10)

MS8132-20 20 sticks (500 mL×20)

Storage: Store at room temperature for 36 months. The product is shipped at room temperature.

Introduction

This product is provided as white or pale yellow instant granules. The main components are Tris, MOPS, EDTA, and SDS. Each stick yields 500 mL of 1× Tris/MOPS/SDS electrophoresis buffer, containing 50 mM Tris, 50 mM MOPS, 1 mM EDTA, and 0.1% SDS.

Quick Start Protocol

Dissolve one stick in distilled water and bring the final volume to **500mL** to obtain 1× Tris/MOPS/SDS electrophoresis buffer.

