



YamayBio

TBS / Tween Buffer Instant Granules

USER GUIDE

For research use only.
Not for use in diagnostic procedures.



Contents and Storage

TBS / Tween Buffer Instant Granules

MS8133 10 sticks (1 L×10)

Storage: Store at room temperature for 36 months. The product is shipped at room temperature.

Introduction

This product is supplied in convenient white or off-white instant granule sticks containing a blend of sodium chloride, potassium chloride, Tris, and Tween-20. Reconstitution of one stick yields 1 liter of 1× TBST buffer with a final concentration of 137 mM sodium chloride, 2.7 mM potassium chloride, 25 mM Tris, and 0.05% Tween-20.

Quick Start Protocol

Dissolve one stick in distilled water and bring the final volume to **1 L** to prepare 1× TBS / Tween Buffer.

