

CERTIFICATE OF ANALYSIS

Product Name: Rapid Coomassie Blue Staining (Destaining-free)
Cat.No: BM3151
Lot.No: A31020
Production Date: 03/05/2025
Expiry Date: 03/05/2026
Storage Temperature: Store at 4°C away from light

Parameter	Specification	Result
Visual appearance	Dark blue, free of particulates	Conforms
Function test	A series of bovine albumin (10ng, 20ng, 40ng, 80ng) were separated by electrophoresis on a 12.5% Tris-glycine SDS -PAGE, 1.0mm, 10wells and then stained with BM3151 Rapid Coomassie Blue Staining Solution	10ng BSA was detected

Quality authorized by:

Esther Zhang

Issue date: 03/05/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com