

CERTIFICATE OF ANALYSIS

Product Name: Cell ExM Kit
Cat.No: CC6232
Lot.No: 64B2110
Production Date: 02/05/2026
Expiry Date: 02/05/2027
Storage Temperature: Store at 2~8℃

Parameter	Specification	Result
Visual appearance	ExM buffer A: colorless, free of particulates	Conforms
	ExM buffer B: light yellow, free of particulates	
	ExM buffer C: colorless, free of particulates	
	ExM buffer D: colorless, free of particulates	
Volume	ExM buffer A : 4-4.5ml/vial	Conforms
	ExM buffer B: 4-4.5ml/vial	
	ExM buffer C: 40-40.5μl/vial	
	ExM buffer D: 40-40.5μl/vial	
Green glass slides	10	Conforms
Glass bottom dishes	4	Conforms
Function test	Hela cells were imaged before and after expansion. ZO-1 (green) and Claudin-1 (red) were detected	Hela cells were expanded 4 times. Colocalization between ZO-1 and Claudin-1 was detected using an OLYMPUS CKX53SF 40x microscope

Quality authorized by: *Esther Zhang*

Issue date: 02/06/2026

YamayBio Inc.
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA
Website: www.yamaybio.com
Email: info@yamaybio.com