

CERTIFICATE OF ANALYSIS

Product Name: Animal Membrane and Cytosol Protein Extraction Kit
Cat.No: BF7475-10
Lot.No: 05B2102
Production Date: 02/03/2026
Expiry Date: 02/03/2027
Storage Temperature: Store Cell Permeabilization Buffer and Solubilization Buffer at -20°C to avoid freeze/thaws.
Store Cell Wash Solution at 4°C.

Parameter	Specification	Result
Visual appearance	Buffer A: colorless, transparent, free of particulates	Conforms
	Buffer B: light yellow, free of particulates	
	Buffer C: colorless, transparent, free of particulates	
Volume	Buffer A: 45-46ml	Conforms
	Buffer B: 10-11ml	
	Buffer C: 5-6ml	
Function test	Membrane proteins and cytosolic proteins were isolated from 5×10^6 HEK293 cells following membrane and cytosol protein extraction kit protocol. GAPDH and Caveolin-1 were analyzed by Western blot.	GAPDH was detected in cytosolic extract while not in membrane extract, and Caveolin-1 was detected in membrane extract while not in cytosolic extract

Quality authorized by:



Issue date: 02/04/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com