

CERTIFICATE OF ANALYSIS

Product Name: Ponceau S Staining Solution
Cat.No: BM3156
Lot.No: 88B2105
Production Date: 02/02/2026
Expiry Date: 02/02/2027
Storage Temperature: Store at 2~8°C

Parameter	Specification	Result
Visual appearance	Purplish-red liquid, free of particulates	Conforms
Function test	For each well of the SDS-PAGE gel, 2 µg, 1 µg, 500 ng, and 200 ng of BSA protein were loaded, along with 3 µl of CY5462 non-prestained protein marker. After electrophoresis, the proteins were transferred to an NC membrane, and all bands were successfully observed following BM3156 staining.	Conforms

Quality authorized by: *Esther Zhang*

Issue date: 02/02/2026

YamayBio Inc.
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA
Website: www.yamaybio.com
Email: info@yamaybio.com