

Protein-Free Blocking Buffer(1x)

Key Features



Rapid and Efficient

Complete membrane blocking in 5 -10 minutes



Protein - Free

Eliminates cross-reactivity with antibodies or target proteins.

High signal-noise ratio provides superior background reduction to non-fat milk and BSA.

Suitable for the detection for phosphorylated proteins.



Cat. No.

Size

BM3152 - 500

500 mL

Applicable to phosphorylated proteins



Higher signal-noise ratio compared to non-fat milk and BSA



Target: HSP60(Cat.No P900009), dilution at 1:5000

MW: 63kD

Lysates: HeLa whole cell lysates

Secondary antibody: Goat anti-Rabbit IgG(H+L),HRP conjugated, diluted 1:5000