

# CERTIFICATE OF ANALYSIS

**Product Name:** ECL Western Blotting Ultimate Sensitivity Substrate (Femto)  
**Cat.No:** VQ4726-100  
**Lot.No:** 54B3110  
**Production Date:** 03/18/2026  
**Expiry Date:** 03/18/2027  
**Storage Temperature:** Stored at 4°C away from light

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
Function test	Serially diluted HRP were dispensed into a 96-well black microplate. After adding VQ4726-100 working solution, distinct concentration-dependent signals were clearly developed within recommended exposure time	Conforms

Quality authorized by: *Esther Zhang*

**Issue date:** 03/20/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)