

CERTIFICATE OF ANALYSIS

Product Name: ECL Western Blotting Sensitivity Substrate (Pico)
Cat.No: VQ4727-500
Lot.No: A2220
Production Date: 11/07/2025
Expiry Date: 11/07/2026
Storage Temperature: Stored at 4°C away from light

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
Function test	Serially diluted HRP were dispensed into a 96-well black microplate. After adding VQ4727-100 working solution, distinct concentration-dependent signals were clearly developed within recommended exposure time	Conforms

Quality authorized by: *Esther Zhang*

Issue date: 11/10/2025

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com