

CERTIFICATE OF ANALYSIS

Product Name: One-Step Stain-free PAGE Gel Rapid Preparation Kit (12.5%)
Cat.No: RF1524
Lot.No: 42B5102
Production Date: 05/06/2026
Expiry Date: 05/06/2027
Storage Temperature: Store at 4°C, the enhanced catalyst is stored at -20°C.

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
pH	8.6-9.0	Conforms
Gelation time	10-30min	Conforms
Migration	the band of 15kd of marker was separated	Conforms
Gel visualization	5µg, 10µg, 20µg and 30µg of total proteins were loaded, then the gels were excited at 302 nm UV light, protein bands can be visualized in 5mins	Conforms

Quality authorized by:

Esther Zhang

Issue date: 05/07/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com