

# CERTIFICATE OF ANALYSIS

**Product Name:** Stain Free Gradient-like PAGE Gel Kit  
**Cat.No:** RF1541  
**Lot.No:** 46B5102  
**Production Date:** 05/06/2026  
**Expiry Date:** 05/06/2027  
**Storage Temperature:** Store at 4°C, the enhanced catalyst is stored at -20°C.

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
pH	8.6-9.0	Conforms
Gelation time	10-30min	Conforms
Migration	the band of 10kd of marker was separated	Conforms
Gel visualization	5µg, 10µg, 20µg and 30µg of total proteins were loaded, then the gels were excited at 302nm UV, protein bands can be visualized	Conforms

Quality authorized by: *Esther Zhang*

**Issue date:** 05/07/2026

YamayBio Inc.  
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA  
Website: [www.yamaybio.com](http://www.yamaybio.com)  
Email: [info@yamaybio.com](mailto:info@yamaybio.com)