

# CERTIFICATE OF ANALYSIS

**Product Name:** One-Step Stain-free PAGE Gel Rapid Preparation Kit (7.5%)  
**Cat.No:** RF1522  
**Lot.No:** 40B5102  
**Production Date:** 05/06/2026  
**Expiry Date:** 05/06/2027  
**Storage Temperature:** Store at 4°C, the enhanced catalyst is stored at -20°C.

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
pH	8.6-9.0	Conforms
Gelation time	10-30min	Conforms
Migration	the band of 35kd of marker was seperated	Conforms
Gel visualization	5µg, 10µg, 20µg and 30µg of total proteins were loaded, then the gels were excited at 302 nm UV light, protein bands can be visualized in 5mins	Conforms

Quality authorized by:

*Esther Zhang*

**Issue date:** 05/07/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)