

CERTIFICATE OF ANALYSIS

Product Name: Cell ExM Plus Kit
Cat.No: CC6233
Lot.No: 65B3102
Production Date: 03/05/2026
Expiry Date: 03/05/2027
Storage Temperature: Stored at 4°C

Parameter	Specification	Result
Visual appearance	ExMPro Reagent A: white powder	Conforms
	Regent A Diluent: light yellow, transparent, free of particulates	
	ExM Reagent B (100×): colorless, transparent, free of particulates	
	ExM Reagent C (100×): colorless, transparent, free of particulates	
Volume	ExMPro Reagent A: 0.51g × 4	Conforms
	Regent A Diluent: 5ml	
	ExM Reagent B (100×): 60μl	
	ExM Reagent C (100×): 60μL	
Casting Frame-C	2 molds: one each of 10mm and 15mm	Conforms
Stickers-C	30 sheets	Conforms
Function test	Hela cells were imaged using Cell ExM Plus Kit before and after expansion. And ZO-1 (green) and Claudin-1 (red) were detected	Hela cells were expanded 4.5 times. Colocalization between ZO-1 and Claudin 1 was detected using an OLYMPUS CKX53SF 40x microscope

Quality authorized by:

Esther Zhang

Issue date: 03/10/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com