

CERTIFICATE OF ANALYSIS

Product Name: ECL Western Blotting Ultimate Sensitivity Substrate (Femto)
Cat.No: VQ4726-100
Lot.No: 54AD1020
Production Date: 08/01/2026
Expiry Date: 08/01/2027
Storage Temperature: Stored at 4°C away from light

Parameter	Specification	Result
Appearance	colorless transparent liquid	Conforms
Function test	Serially diluted HRP were dispensed into a 96-well black microplate. After adding VQ4726-100 working solution, distinct concentration-dependent signals were clearly developed within recommended exposure time	Conforms

Quality authorized by:

Xue Zhu

Issue date: 08/03/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: www.yamaybio.com

Email: info@yamaybio.com