

CERTIFICATE OF ANALYSIS

Product Name: Tricolor Prestained Protein Ladder 10 kDa~250 kDa (Orange, Green & Blue)
Cat.No: CY5466-2500
Lot.No: 18B71020
Production Date: 07/06/2026
Expiry Date: 07/06/2028
Storage Temperature: Store at 4°C for 3 months or store at -20°C for 24 months

Parameter	Specification	Result
Band Pattern and mobility	Separation by 15% SDS-PAGE and post-electrophoretic in-gel analysis of band migration. Transferred on the membrane after electrophoresis	11 bands were sharp and flat. The migration of the 11 bands met the product specifications

Quality authorized by: *Esther Zhang*

Issue date: 07/07/2026

YamayBio Inc.
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA
Website: www.yamaybio.com
Email: info@yamaybio.com