

# CERTIFICATE OF ANALYSIS

**Product Name:** Tricolor Prestained Protein Ladder 10 kDa-250 kDa (Orange, Green & Blue)  
**Cat.No:** CY5466-500  
**Lot.No:** 18B71002  
**Production Date:** 06/01/2026  
**Expiry Date:** 06/01/2028  
**Storage Temperature:** Store at 4°C for 3 months or store at -20°C for 24 months

| Parameter                 | Specification                                                                                                                               | Result                                                                                     |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
| Band Pattern and mobility | Separation by 15% SDS-PAGE and post-electrophoretic in-gel analysis of band migration.<br>Transferred on the membrane after electrophoresis | 10 bands were sharp and flat. The migration of the 10 bands met the product specifications |

Quality authorized by:

Xue Zhu

**Issue date:** 06/02/2026

YamayBio Inc.

Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA

Website: [www.yamaybio.com](http://www.yamaybio.com)

Email: [info@yamaybio.com](mailto:info@yamaybio.com)