

# CERTIFICATE OF ANALYSIS

**Product Name:** Tris/Glycine/SDS Buffer Instant Granules  
**Cat.No:** MS8131-10  
**Lot.No:** 47B71100  
**Production Date:** 07/2026  
**Expiry Date:** 07/2029  
**Storage Temperature:** Stored at room temperature

| Parameter         | Specification                                       | Result   |
|-------------------|-----------------------------------------------------|----------|
| Visual appearance | White or light yellow granules                      | Conforms |
| Solubility        | 200s(500rpm)                                        | Conforms |
| Weight            | 17.52-19.36/pack                                    | Conforms |
| pH                | 8.1-8.5                                             | Conforms |
| Conductivity      | 600-800 $\mu$ S/cm                                  | Conforms |
| Dnase             | none                                                | Conforms |
| Protease          | none                                                | Conforms |
| Function test     | The detected bands have high resolution and clarity | Conforms |

Quality authorized by: *Xue Zhu*

**Issue date:** 07/2026

YamayBio Inc.  
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA  
Website: [www.yamaybio.com](http://www.yamaybio.com)  
Email: [info@yamaybio.com](mailto:info@yamaybio.com)