CERTIFICATE OF ANALYSIS

Product Name:	PBS Buffer Instant Granules
Cat.No:	MS8138
Lot.No:	74B71100
Production Date:	07/2026
Expiry Date:	07/2029
Storage Temperature:	Stored at room temperature

Parameter	Specification	Result
Visual appearance	white particles	Conforms
Solubility	100 s (500 rpm)	Conforms
Weight	9.11-10.07 g/pouch	Conforms
PH (1 pouch/1L Water）@25℃	7.2-7.6	Conforms
Conductivity (1 pouch/1L Water）@25℃	13.0-16.5 mS/cm	Conforms

Quality authorized by: Xue Zhu

Issue date: 07/2026

YamayBio Inc.
Address: 6191 Cornerstone Court, Suite 107, San Diego, CA 92121, USA
Website: www.yamaybio.com
Email: info@yamaybio.com